IDitalian *design*

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EN PRECAUTIONS

Please read the instructions carefully before using the appliance. To reduce any risk of accident, take the following precautions:

1. Always unplug the appliance immediately after use.
2. Never use the appliance while you are in the bath.
3. IEC 60417 - 5582: Do not place or store the appliance anywhere where it could fall or be pushed towards the bath or sink. 
4. Do not immerse or drop into the water or any other liquid.
5. If the appliance falls into the water, unplug it immediately. Do not place your hand in the water.
6. This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
7. Keep the appliance away from water.
8. An appliance should never be left unattended while it is plugged into the mains.
9. Use this appliance only for its intended purpose, as indicated in this manual. Do not use any accessory that is not recommended by the manufacturer.
10. Do not use the appliance if you are drowsy.
11. Do not drop the appliance or insert any object in the openings or tube.
12. If your skin is under nursing, or under any unhealthy condition, please consult your doctor before using the device.

13. Clean the residue of gels or creams with a soft damp cloth or paper towowel immediately after each use.

14. This care instruments suitable for private use only, do not share with your family or friends.

15. Preserve the device in a cool, dry place.

16. Before cleaning the instrument, make sure the power is off.

17. This product must not be used by children or by other people without assistance or supervision if their physical, sensory or mental condition limits safe use of the appliance.

Children must be watched to make sure that they do not play with the appliance.

18. Do not use the appliance in case of failure and/or malfunction. Contact directly the technical service.

EN PROPERTIES

Ultrasonic wave is a kind of high frequency mechanical vibration wave, which is widely used in medical and beauty treatments fields. The mechanical wave has mechanical, thermal and catalytic effects. It can produce mild resonance effect on human skins when it works on low energy and high frequency vibration, enhance cell vitality and permeability between cortical tissues, dredge sweat glands and contract pores.

The ultrasonic deaging and beauty instrument uses the ultrasonic and nano technology, make an intense friction between molecules with a vibration of 24000 per second, and the vibration can penetrate into the skin to make the intense friction between the molecules, the Deep skin temperature is raised, the metabolism of the cells is strong, the microcirculation is accelerated, and the return of the metabolic products is promoted. Through the transformation of electron flow, the permeability of cell membrane is enhanced, which make the semi-permeable membrane produce intracellular primary plasma microflow, and the volume of fatty acid alienation into CO₂ and H₂O fat is reduced. At the same time, with a certain target product, the beauty effect is more remarkable.

EN FUNCTIONS

1. Deep cleaning, clear aged cutin.
2. Decompose pigment and reduce the color spots
3. Remove wrinkle in Deep, increase skin elasticity.
4. Promote blood circulation and lymphoid metabolism.
5. Enhance the regeneration ability and permeability of skin cells.
6. Massage skin, nutrition leading in and mixed color leading-out.
7. No damage and side effects to the skin
8. Beautiful appearance, easy to use, simple operation and easy to carry.

EN PRODUCT PROBLEMS

Symptom	Reason	Processing Method
The product can't be turned on	Undercharge	Please refer to the instruction charging
	Less of power electricity	
Shut down in use	Automatic power off function (automatic cut off function about 20 minutes)	It is normal function, not caused by failure
	Less of power electricity (the product will beep and beep during operation)	Please refer to the usage of the instructions
Unable to charge	The USB interface is not connected	Please plug in the USB cable
	The lithium battery charging device for mobile phone used in charging	Please replace the charging device
Unable to charge	The charging device has a fault	Please restart the charging device
	Adapter (sold separately) failure	Please use another charging device

EN CHARGING METHOD

1. USB Charging line connect to USB interface. The battery is not full charged before shipment, please charge it before use.
2. Plug in the USB interface of a computer, an adapter (not included), or a lithium battery of charge of a mobile phone. The LED light is continuous flashing, and the LED light is on after the charging is finished, please confirm the status of the lithium battery before charging.

EN CHARGING ATTENTIONS

Please charge it if the product dose not display the boot status when long press the power button.
The product will make a "beep-beep" warning sound in operation, then the shutdown state will appear. It shows that the battery's margin is insufficient, please charge the product immediately. The product is automatically standby after working 10 minutes, and if it is not used for one minute, it will be turned off automatically. It is about 120 minutes to be fully charged. (The first use requires charging more than 5 hours)
It can be used for about 60 minutes after full charged.
It is recommended that the product should be periodically charged.

EN OPERATION GUIDE

Cleansing mode: Long press the power button 1.5 seconds to turn on the machine.

After the machine is powered on, it will automatically start the first cleansing mode. Daily cleansing: Black head, Acne, Cutin. Use with facial cleanser, ultrasonic penetration onto muscle floor, deeply cleansing pores.

Moisturizing mode: Long press the power button 1.5 seconds to turn on the machine.

Short press the power button one time to start the second moisturizing mode. The ion leading-in the nutrition, the essence is input into the bottom of the muscle under the action of the ultrasonic wave and the negative ion, and the absorption is accelerated.

Note: The switch lights will flash when the mode changed.

Lifting mode: Long press the power button 1.5 seconds to turn on the machine

Short press the power button 2 times to start the third lifting mode EMS microcurrent, restore the elasticity of the skin by muscle massage through microcurrent stimulation. At the same time, EMS and the ultrasonic wave transmits the microcurrent to the muscle bottom through the inductor, the collagen is regenerated, so that the effect of the compact lifting and the beauty treatments are archived.

EN MAINTENANCE AND STORAGE

1. Make sure the machine is turned off.
2. Please wipe off the stains on the sholvel head and the unit body by a wet paper towel or alcohol. Do not wipe the product with corrosive chemical industrial products such as diluent, solvent, etc.
3. Please put it in a dry place for storage after cleaning. Do not keep it in wet places for storage, such as bathrooms, etc. . . Do not keep it in a place where the children can touch.

EN ENVIRONMENT

Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries with separate collection programs. This symbol, found on your product or on its packaging, indicates that this product should not be treated as

household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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